



## Our new wheelchair mounting solution for your arm support

Rehadapt's tireless commitment to expand boundaries of interaction - literally in every direction - led us to a new mounting solution for the Edero arm support system from Armon Products. It's not just about supporting an arm. Nor is it just about mounting another product. It's about reaching a further step of independence.



The **Armon Edero** is an arm support, especially designed for those with limited strength in their arms, or for those that can only use their muscular strength for a limited time.

The device compensates the full weight of the arm with an innovative spring compensation mechanism, easy to setup and intuitive to use. The Edero follows the natural movement of the user and allows the arm to move easily in all directions.

That gives the beneficiary full control over his arm movements.

## M3D Edero Arm Support Adapter 16.1022



Rehadapt's device adapter for mounting the Edero Arm Support from Armon includes a stainless steel swivel interface to carry the Edero and a hollow element to park it next to the armrest. Both elements are ready to fit perfectly into our M3D ecosystem of mounting elements.

With this adapter, the Edero Arm Support can be mounted onto any wheelchair:



under the armrest of a  
powered wheelchair



on the slide-track of a  
powered wheelchair



anywhere on the frame  
of a manual wheelchair



With our unique virtual mounting service it will be even easier for you to find the right mounting solution for your Edero arm system. This option will be integrated into our VMS process soon.



For more information feel free to contact our [salesteam@rehadapt.com](mailto:salesteam@rehadapt.com)  
or visit [www.rehadapt.com/armon](http://www.rehadapt.com/armon)