




# Sockets

"The mother of all sockets" - our universal device socket (UDS) includes our quick release mechanism which allows for the safe device attachment within seconds. The footprint on the device is minimal, yet great sturdiness is achieved!



# Sockets



## UDS 22 DC | 16.1051

Rehadapt's signature quick release system with ball-joint in three versions.

- ✓ with 3 axes of rotation
- ✓ Two levers for separated fixation of the ball joint and the tube clamp. This way, the weight of the device remains supported during repositioning.
- ✓ The included bolts and Quickshifter-lever allow you to replace or optimize the existing fastening elements.




## UGA 22 | 16.1000

- ✓ Fixation with a bolt.



## UGA 22 QS | 16.1050


- ✓ Fixation with a lever.



## UDS Eco 16.1004

- ✓ with 1 axis of rotation

The most economic and compact version of Rehadapt's quick release system: sitting directly on the tube. Allows for the rotation around the tube axis. Fastened with a bolt.



## UDS Eco Plus 16.1005

- ✓ with 2 axes of rotation

Excellent compromise between economic yet flexible use of Rehadapt's quick release system. Allows for the rotation around the tube axis and the twisting of the screen. Fastened with a bolt.



## System bridges

Combine your existing mounting elements with Rehadapt's ecosystem.



**SLS interface**  
→ to UDS  
14.1516



**VESA fitting**  
→ to UDS  
16.1010

Perfect for 3rd party monitor holders to carry devices with a Rehadapt UDS adapter plate.



**UDS interface**  
→ to Daessy  
16.1020



**EyeGaze  
Edge Talker**  
16.1021



**SLS interface**  
→ to 1/4"-UNC  
14.9511\_SW



**Edero arm support**  
→ to M3D  
16.1022

Device adapter for mounting the Edero Arm Support from Armon on Rehadapt's Monty 3D system. Includes a pin for inserting the arm and a tube for parking it.



**Universal  
70/50**  
16.1224



## Wheelchair mounting solution for your arm support


Rehadapt's tireless commitment to expand boundaries of interaction – literally in every direction – led us to a new mounting solution for the Edero arm support system from Armon Products. It's not just about supporting an arm. Nor is it just about mounting another product. It's about reaching a further step of independence.



The Armon Edero is an arm support, especially designed for those with limited strength in their arms, or for those that can only use their muscular strength for a limited time.

The device compensates the full weight of the arm with an innovative spring compensation mechanism, easy to setup and intuitive to use. The Edero follows the natural movement of the user and allows the arm to move easily in all directions.

That gives the beneficiary full control over his arm movements.



## M3D Edero Arm Support Adapter 16.1022

Rehadapt's device adapter for mounting the Edero Arm Support from Armon includes a stainless steel swivel interface to carry the Edero and a hollow element to park it next to the armrest. Both elements are ready to fit perfectly into our M3D ecosystem of mounting elements.

With this adapter, the Edero Arm Support can be mounted onto any wheelchair:



under the armrest of a  
powered wheelchair




on the slide-track of a  
powered wheelchair



anywhere on the frame of  
a manual wheelchair



With our unique virtual mounting service it will be even easier for you to find the right mounting solution for your Edero arm system. This option is integral part of our VMS process. Learn more at [www.rehadapt.com/vms](http://www.rehadapt.com/vms)



For more information feel free to contact our [salesteam@rehadapt.com](mailto:salesteam@rehadapt.com)  
or visit [www.rehadapt.com/armon](http://www.rehadapt.com/armon)