


Technical Reference

Floorstands



Rehadapt Engineering
GmbH & Co. KG

📖 Heinrich-Hertz-Str. 104
34123 Kassel, Germany

☎ Tel.: +49 561 220 717 0

📠 Fax: +49 561 220 717 99

✉ officeteam@rehadapt.com

🌐 <https://rehadapt.com>

Rehadapt
North America

📖 40 Oser Ave, Suite 40-01
Hauppauge, NY 11788


☎ Tel.: 1 (888) 635-6962

📠 Fax: 1 (888) 635-6962


✉ officeteam.na@rehadapt.com

🌐 <https://rehadapt.com>

Content





17.1055	FS EcoLock	3
17.1050	FS TeleLock	4
17.1052	FS TeleLock QP	5
17.1061	FS VarioLock	6
17.1063	FS VarioLock QP	7
17.1093	FS VarioLock QP Plus	8
17.1080	FS EcoFloat Light	9
17.1070	FS EcoFloat	9
17.1072	FS EcoFloat HD	9
17.1073	FS TeleFloat	10
17.1083	FS TeleFloat Light	10
17.1074	FS TeleFloat QP	11
17.1084	FS TeleFloat Light QP	11
17.1065	FS VarioFloat	12
17.1085	FS VarioFloat Light	12
17.1066	FS VarioFloat QP	13
17.1086	FS VarioFloat Light QP	13







17.1102	FS PentaLock Eco	14
17.1100	FS PentaLock Tele	15
17.1172	FS PentaFloat Eco	16
17.1170	FS PentaFloat Tele	17







17.1400	FS Mini	18
17.1401	FS Mini SLS	18
17.1410	FS Mini HD	19
17.1411	FS Mini SLS HD	19














370 mm
[14.5 inch]

400 mm
[16 inch]

FS VarioLock QP Plus

17.1093




1130 - 1830 mm
[44.5 - 72 inch]

9.3 kg
[20.6 lbs]

77 mm
[3 inch]

670 - 1175 mm
[26.4 - 46.3 inch]

780 mm
[30.7 inch]





FS EcoFloat Light


17.1080

FS EcoFloat

17.1070

FS EcoFloat HD


17.1072




	17.1080	17.1070	17.1072*
	0.7-3.0 kg [1.5-7 lbs]	2-6 kg [4.4-13.2 lbs]	6-11 kg [13.2-24.2 lbs]




1x
17.9610
17.9611
17.9614



1x
17.9608



1x
17.9530 / 17.9533



2x
17.9500

4x
01.0850

1x
01.0620

1x
01.0657

1x
01.0436

1x
01.0437

1x
01.0027





FS TeleFloat


17.1073

FS TeleFloat Light


17.1083




1x
17.9610
17.9614



1x
17.9609



1x
17.9530




2x
17.9500



4x
01.0850



1x
01.0620



1x
01.0657



1x
01.0436



1x
01.0437





1x
01.0027

FS TeleFloat QP

17.1074

FS TeleFloat Light QP



17.1084





FS VarioFloat 17.1065


FS VarioFloat Light 17.1085







FS VarioFloat QP 17.1066

FS VarioFloat Light QP 17.1086










	17.1086	17.1066
 max.	0,7-3,0 kg [1.5-7 lbs]	2-6 kg [4.4-13.2 lbs]













FS PentaLock Eco

17.1102





FS PentaLock Tele 17.1100







FS PentaFloat Tele

17.1170



FS Mini
17.1400

FS Mini SLS
17.1401





* 17.1400

FS Mini HD

17.1410

FS Mini SLS HD

17.1411



* 17.1410

