

Bases

11.1060 - M3D RKL E30-70



2024-02-01
00.5000-H

Rehadapt GmbH


Heinrich-Hertz-Str. 104
34123 Kassel, Germany
Tel.: +49 561 220 717 0
Fax: +49 561 220 717 99
officeteam@rehadapt.com
<https://rehadapt.com>

MADE EASY MOUNTING




Rehadapt North America


40 Oser Ave, Suite 40-01
Hauppauge, NY 11788
Tel.: 1 (888) 635-6962
Fax: 1 (888) 635-6962
officeteam.na@rehadapt.com
<https://rehadapt.com>




- 4x 01.0626
● M6x25
- 4x 01.0640
● M6x40
- 4x 01.0614
● M6x55



X	
0-5 mm [0-0.2"]	01.0626
5-20 mm [0.2-0.8"]	01.0640
21-40 mm [0.8-1.55"]	01.0614



CAUTION / VORSICHT



NOTE / HINWEIS

