

Adapters

14.1642 -
L3D GA Q-Logic 3/ LiNX 500/ Omni 2

00.5282-A 2024-11-15



Rehadapt GmbH


Heinrich-Hertz-Str. 104
34123 Kassel, Germany
Tel.: +49 561 220 717 0
Fax: +49 561 220 717 99
officeteam@rehadapt.com
https://rehadapt.com

MADE EASY



Rehadapt North America

40 Oser Ave, Suite 40-01
Hauppauge, NY 11788
Tel.: (888) 635-6962
Fax: (888) 635-6962
officeteam.na@rehadapt.com
https://rehadapt.com




1x

2x 01.0566
● M4x10



1x 01.0002
● 3 mm

3x 01.0157
● M5x14





1x 01.0035
● 2.5 mm



44 mm
[1.73"]

30 mm
[1.18"]




105 mm
[4.1"]




Mounting options / Montagemöglichkeiten

Option 1

(LiNX 500)




2x 01.0157
M5x14




Option 2

(Q-Logic 3 EX Enhanced Display)




3x 01.0157
M5x14




Option 3


(R-net Omni 2)




2x 01.0566
M4x10



Option 1
(LiNX 500)



Option 2
(Q-Logic 3 EX Enhanced Display)



Option 3
(R-net Omni 2)

