

Tubes

M3D



M3D Tubekit
#11.1400

Monty-3D Standard tube kit, including 3 tubes in lengths of 200mm (8"), 400 (16") and 500mm (20").



R22,2-200 to R22,2-500

Tube, 22,2 mm diameter with following lengths:

- 200 mm (8") #01.5005
- 300 mm (12") #01.5010
- 400 mm (16") #01.5020
- 500 mm (20") #01.5030



Slip-Out-Safeguard

R22,2-200 SOS to R22,2-500 SOS

Tube with Slip-Out Safeguard, 22,2 mm diameter with following lengths:

- 200 mm (8") #01.5505
- 300 mm (12") #01.5510
- 400 mm (16") #01.5520
- 500 mm (20") #01.5530




Slip-Out Safeguard kit
#01.5500

2 slip-out safeguards for tubes Ø22x2.3mm with screws. Prevents tubes from slipping out. (without allen wrench)

S-Tubes



S-tubes in different sizes with offset for extra clearance.



M3D-HD



M3D-HD Upgrade-kit
#11.1403

Monty-3D upgrade kit to equip any M3D mount with a double tube. Significantly reduces vibration, allows for heavy devices on long mounts.



Easy-Off Kit for EyeControl HD
#11.1510

Allows for quick removal of the entire mount directly at the wheelchair frame. (Only for double tube mounts)



M3D-Connector HD
#11.1500

Monty-3D upgrade kit to equip any M3D mount with a double tube. Significantly reduces vibration, allows for heavy devices on long mounts.

Tubes

L3D



**R10-100
to R10-350**

**R10-100
to R10-350**



R10-50
#14.9511



L3D-Rod, 10mm diameter with following lengths:

L3D-Rod, 16mm diameter with following lengths:

L3D-Rod, Ø 10 x 50 mm (3/8" x 2") with 1/4" threading for switch adapters on one side.

100 mm (4")	#14.9512
150 mm (6")	#14.9513
200 mm (8")	#14.9514
250 mm (10")	#14.9515
300 mm (12")	#14.9516
350 mm (14")	#14.9517

100 mm (4")	#14.9500
150 mm (6")	#14.9501
200 mm (8")	#14.9502
250 mm (10")	#14.9503
300 mm (12")	#14.9504
350 mm (14")	#14.9505
400 mm (16")	#14.9506

