

# How to order a Light-3D Mount?

## A Where to mount to?



Table or Board



Wheelchair



Power Wheelchair

Select products with element „TC“ in product name or product numbers starting with #14.4xxx.

Select products with element „WC“ in product name or product numbers starting with #14.3xxx. Comes with 25/22mm (1”/ 7/8”) round frame clamp by default, please specify any other needed diameter on order.

If other than round frameclamp is required for the chair, please specify on order. Custom frameclamp is included free of charge!

Example: L3D-TC 1AK

Example: L3D-WC 1AK



TC= Table clamp  
WC= Wheelchair

## B How many tubes are needed?



One tube

Product name contains number „1“ after TC / WC.

Example:  
L3D-WC 1AK



Two tubes

Product name contains number „2“ after TC / WC.

Example:  
L3D-WC 2AK



Three tubes

Three-tube-System always includes a swivel joint. Product name contains number „3“ after TC / WC.

Example:  
L3D-WC 3AK

## C Joints with Screws or Clamping Levers?



Allen Key Screw

Product name contains element „AK“ after the number.

Example:  
L3D-WC 1AK



Quick Shift Lever

Product name contains element „QS“ after the number.

Example:  
L3D-WC 2QS



AK= Allen Key  
QS= Quick Shift

## D Pin or quick release connector?



SLS

By default Light-3D mounts end with a joint to chuck a pin from the device cradle („SLS= Spigot Link System”).

Example: L3D-WC 2QS



UGA

If a quick release connector is desired („UDS= Universal Device Socket”), please order products containing „U“ in the name.


Example: L3D-WC 1AK-U

## E Order Device Holder SEPARATELY



SLS

A device holder needs to be ordered **separately**. According to prior decision choose from Article Numbers #14.15xx (with SLS-Spigot) or #16.11xx (with UDS connector).



UGA